Ionplus⁵⁵

Task List

2083 - Single Run

2072 - Auto Tuning - ANA Steerer Y

2073 - Auto Tuning - LE Magnet

✓ Successful 2021-11-12 13:27

✓ Successful 2021-11-12 13:27

Successful 2021-11-12 14:42

2075 - Auto Measurement - DEMO58

Successful 2021-11-12 16:05

2076 - Sleep

Successful 2021-11-12 16:06

2084 - Wake Up - Immediately

△ COLLAPSE HISTORY

Cycle time 17 of 20 s

Run AC1636731346

Target code 0x20FEF0C3

1 of 12

2085 - Auto Measurement - DEMO58

3 of 5

Sequence 2 of 5

Target 42.4711.114

ılı.

莊

.

ACS USER INTERFACE

SELECT YOUR VIEW

Each view contains a number of cards with specific information and controls to help you run your tuning, visualize data, do maintenance or change settings.

TASK LIST

Schedule tasks such as wake up, tuning, measurement and sleep. Let a wizard task guide you through maintenance. A task planner allows you to save your favorite autotuning and measurement routines.

TUNING CARD

Access all devices relevant to the ion beam. Tune manually or select auto-tuning for the fastest tuning of a device. Save tuning sets and load settings from previous measurements.

27767

-1000

24702

Stripper Contro

HE Magnet


HE ESA

Tuning Settings

CHAT & LOG, NOTIFICATIONS

Get an overview of current and past events. Comment on tasks and notify other users.

Measuring - Position: 3 - Sequence: 2 Vacuum




-0.31 A

37.69 kV -37.69 kV

Q Filter by tag

-0.10 V

Sample Changer Camera



PLOT CARD

Visualize data and look at tunings There is also a corresponding view that allows you to look at graphs and tables of all runs.

MANAGEMENT

USER

Each user has their own account and only one user at a time is 'in control'.

STATUS CARDS

See the most important instrument information such as currents, pressures and device states at a glance.